

## INSTRUCTION FOR USE

### Preva EDTA (17% EDTA Solution)

Preva EDTA is a 17% Ethylenediaminetetraacetic Acid (EDTA) solution buffered to neutral pH, intended for use during root canal treatment. It facilitates rapid cleansing, preparation, and enlargement of root canals through chelating action and assists in smear layer removal during endodontic procedures.

**COMPOSITION**  
Disodium Edetate (17%), Purified Water, Excipients

**INTENDED USE**  
Preva EDTA is intended for use as an irrigating solution during root canal treatment procedures.

- INDICATIONS FOR USE**
- Root canal cleansing during endodontic procedures
  - Root canal shaping during instrumentation
  - Removal of smear layer from root canal walls
  - Facilitation of opening and negotiation of calcified canals
  - Chelation and softening of dentin during root canal preparation

- CONTRAINDICATIONS**
- Do not use in patients with known hypersensitivity or severe allergic reactions to any components of the product
  - Do not use for applications other than intended use

- WARNINGS**
- Avoid prolonged or repeated contact with skin and eyes
  - Avoid prolonged inhalation of vapors or aerosols
  - Do not ingest the product
  - Contact with skin or mucosal tissues may cause irritation or hypersensitivity reactions
  - Avoid extrusion of the solution beyond the root apex
  - Keep out of reach of children

In case of accidental exposure:  
Skin contact:

- Immediately remove material using cotton and rinse thoroughly with soap and water.

- Eye contact:
- Immediately flush eyes with copious amounts of clean water and seek medical attention.

- Ingestion:
- Seek immediate medical advice.
- If skin rash, irritation, sensitization, or any allergic reaction occurs, discontinue use and seek medical attention.

- PRECAUTIONS**
- Product intended for professional dental use only
  - Wear appropriate personal protective equipment including gloves, masks, and protective eyewear during handling
  - Use rubber dam isolation whenever possible
  - Avoid injecting the solution beyond the apex of the canal
  - Use only the amount necessary for the intended clinical procedure
  - Do not return dispensed material to the original container
  - Do not reuse dispensed solution after application

#### ADVERSE REACTIONS

- Possible adverse reactions may include:
- Local tissue irritation
  - Hypersensitivity reactions
  - Temporary discomfort upon accidental tissue contact

#### DIRECTIONS FOR USE

1. Transfer the required amount of Preva EDTA solution from the bottle into a disposable syringe fitted with a disposable irrigation needle.
2. Insert the irrigation needle gently into the prepared root canal until slight resistance is felt.
3. Withdraw the needle approximately 1–2 mm to prevent binding within the canal.
4. Slowly irrigate the canal with the solution to facilitate debris removal and dentin chelation.
5. Continue irrigation according to standard endodontic procedures and clinical judgment.
6. Dispose of any remaining dispensed solution after use.

#### Important:

- Do not force the irrigation needle into the canal.
- Do not apply excessive pressure during irrigation.
- Avoid extrusion beyond the apical foramen.

#### PRESENTATION

Ref.No: PL-0183- Preva EDTA – 17% EDTA Solution- 15ml Bottle

#### STORAGE CONDITIONS

- Store at temperatures between 10°C and 24°C (50°F–75°F)
- Keep the container tightly closed after use
- Store in a cool and dry place away from direct sunlight and moisture
- Shelf life is not affected after opening when stored according to instructions
- Do not use after the expiry date


#### SHELF LIFE

3 years from the date of manufacture

#### DISPOSAL














Dispose of unused product, packaging, and contaminated materials according to applicable local regulations and biomedical waste disposal requirements.

#### CONTACT

 **Prevalent Healthcare LLP**, Poongode Road, Karutheni, Malappuram, Kerala – 679339, India.  
Ph.: +91 4931257954 Email: info@prevalentonline.com  
www.prevalentonline.com

Failure to comply with the specified storage conditions may alter the working characteristics of the material and reduce its shelf life. The manufacturer is not responsible for any loss of quality resulting from non-compliance with the specified conditions for transportation, storage, and use of this product. The user assumes full responsibility for using the material for any purposes other than those specified by the manufacturer. This material is intended solely for professional dental use. Application must be carried out strictly according to the instructions for use. The manufacturer accepts no liability for damages resulting from failure to follow the provided instructions or for use beyond the stipulated applications. The user is responsible for assessing the material's suitability for any purpose not explicitly stated in this instruction sheet. The descriptions and data provided do not constitute a warranty of attributes and are not legally binding.

**This product is specially formulated for use in dentistry and to be used by qualified dental personnel**

-  Symbol for "TEMPERATURE LIMITATION"  Symbol for "KEEP DRY"
-  Symbol for "KEEP OUT OF THE REACH OF CHILDREN"
-  Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"
-  Symbol for "KEEP AWAY FROM SUNLIGHT"  Symbol for "CAUTION"
-  Symbol for "MANUFACTURER NAME & ADDRESS"
-  Symbol for "CONSULT INSTRUCTION FOR USE"
-  Symbol for "CATALOGUE NUMBER"  Symbol for "USE BY"
-  Symbol for "MEDICAL DEVICE"
-  Symbol for "NON STERILE"
-  Symbol for "BATCH CODE"

Document No	PL-IFU-029
Revision	02
Date	19-11-2025