



Back to
Home

INSTRUCTION FOR USE

HEMOSTATIC SOLUTION

(25% Buffered Aluminium Chloride Solution)

Prevalent Hemostatic Solution is a 25% Buffered Aluminium Chloride Solution intended for topical use during dental procedures. It is designed to temporarily restrict and control minor gingival bleeding and gingival crevicular fluid in the sulcus surrounding the tooth, helping to provide a clean and dry operating field.

COMPOSITION

Aluminium Chloride – 25%, Glycerine, Purified Water, Excipients

INTENDED USE

Prevalent Hemostatic Solution is intended for use as a topical hemostatic agent to temporarily control minor gingival bleeding and gingival crevicular fluid during dental procedures.

INDICATIONS

Prevalent Hemostatic Solution is indicated for the temporary control of minor gingival bleeding and fluids during:

- Restorative procedures.
- Crown and bridge preparations.
- Impression-taking procedures.
- Oral prophylaxis.
- Other dental procedures where temporary control of minor gingival bleeding or crevicular fluid is required.

WARNINGS

- For professional dental use only.
- The product should be used only by, or under the supervision of, a qualified dental professional.
- For topical use only. Do not inject.
- Do not swallow or take internally.
- Avoid prolonged or repeated contact with skin, eyes, and oral soft tissues.
- Contact with the product may cause irritation or hypersensitivity in susceptible individuals.
- Do not use in patients with a known allergy or hypersensitivity to Aluminium Chloride or any other component of the product.
- Keep out of reach of children.

PRECAUTIONS

- Avoid unnecessary contact with skin, eyes, and mucous membranes.
- Avoid inhalation of mist or vapour generated during handling.
- Do not ingest.
- Avoid excessive or prolonged application to gingival tissues.
- Thoroughly rinse the treated area with water after use.
- Prevent contamination of the contents during dispensing.
- Close the bottle securely immediately after use.



Eye Contact

- The product may cause eye irritation. In case of accidental eye contact, immediately rinse thoroughly with plenty of clean water. Seek medical attention if irritation persists or if clinically indicated.

Skin Contact

- The product may cause skin irritation or sensitization in susceptible individuals. In case of accidental skin contact, immediately remove the material and wash the affected area thoroughly with soap and water.
- If irritation, rash, sensitization, or another allergic reaction develops, discontinue exposure and seek appropriate medical attention.

Accidental Ingestion

- If accidentally swallowed, rinse the mouth thoroughly with water. Do not induce vomiting. Seek medical advice if a significant quantity has been swallowed or if adverse symptoms occur.

CONTRAINDICATIONS

Prevalent Hemostatic Solution should not be used:

- In patients with a known allergy or hypersensitivity to Aluminium Chloride or any other component of the product.
- Where the use of a topical gingival hemostatic agent is clinically inappropriate.
- For the management of significant or uncontrolled surgical bleeding.

INSTRUCTIONS FOR USE

- Isolate the treatment area as clinically appropriate.
- Dispense a small quantity of Hemostatic Solution into a clean disposable container or onto a suitable dispensing surface.
- Saturate a clean cotton pellet or suitable applicator with the solution.
- Apply the saturated cotton pellet directly to the area of minor gingival bleeding requiring hemostasis.
- Maintain contact with the treatment area until adequate hemostasis is achieved.
- Remove the cotton pellet and thoroughly rinse the treated area with water.
- Gently dry and inspect the area to confirm that bleeding and fluid are adequately controlled.
- Proceed with the intended restorative, prosthodontic, impression-taking, or other dental procedure.
- If adequate hemostasis is not achieved, the application may be repeated as clinically necessary. Remove the previous application and rinse the area before reapplication.

Important Application Information

- Use only the minimum quantity required to achieve the desired hemostatic effect.



- Do not allow the bottle or dispensing tip to come into direct contact with oral tissues, contaminated instruments, cotton pellets, or other potentially contaminated surfaces.
- Do not return dispensed material to the original container.
- Material that has already been dispensed or used must not be reused.
- Thoroughly rinse residual solution from the treatment area before proceeding with impression-taking, bonding, or restorative procedures.

STORAGE

- Store at 10°C to 24°C (50°F to 75°F).
- Store in a cool, dry place.
- Protect from excessive heat and direct sunlight.
- Keep the bottle tightly closed when not in use.
- Prevent contamination of the contents.
- Do not use after the expiry date stated on the packaging.

PRESENTATION

Ref. No.: PL-0118- Prevalent Hemostatic Solution – 15 mL Bottle

SHELF LIFE

3 years from the date of manufacture, when stored under the recommended storage conditions.

DISPOSAL

Dispose of unused or expired material, contaminated applicators, cotton pellets, empty containers, and packaging in accordance with applicable local regulations and the prevailing waste-disposal requirements of the country in which the product is used.

CONTACT

Prevalent Healthcare LLP, Poongode Road, Karutheni, Malappuram, Kerala – 679339, India.
Ph.: +91 4931257954 Email: info@prevalentonline.com
www.prevalentonline.com

Failure to comply with the specified storage conditions may alter the working characteristics of the material and reduce its shelf life.

The manufacturer is not responsible for any loss of quality resulting from non-compliance with the specified conditions for transportation, storage, and use of this product. The user assumes full responsibility for using the material for any purposes other than those specified by the manufacturer. This material is intended solely for professional dental use. Application must be carried out strictly according to the instructions for use. The manufacturer accepts no liability for damages resulting from failure to follow the provided instructions or for use beyond the stipulated applications. The user is responsible for assessing the material's suitability for any purpose not explicitly stated in this instruction sheet. The descriptions and data provided do not constitute a warranty of attributes and are not legally binding.



This product is specially formulated for use in dentistry and to be used by qualified dental personnel

Symbol for "TEMPERATURE LIMITATION"

Symbol for "KEEP OUT OF THE REACH OF CHILDREN"

Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"

Symbol for "KEEP AWAY FROM SUNLIGHT"

Symbol for "MANUFACTURER NAME & ADDRESS"

Symbol for "CONSULT INSTRUCTION FOR USE"

Symbol for "CATALOGUE NUMBER"

Symbol for "MEDICAL DEVICE" Symbol for "KEEP DRY"

Symbol for "NON STERILE" Symbol for "USE BY"

Symbol for "BATCH CODE" Symbol for "CAUTION"

Document No	PL-IFU-037
Revision	02
Date	12-01-2026