



Back to
Home

INSTRUCTION FOR USE

Preva Block

(25% Aluminium Chloride Hemostatic Gel)

Preva Block is a 25% Aluminium Chloride Hemostatic Gel intended for topical application in dental procedures. It helps restrict and control minor gingival bleeding and gingival crevicular fluid in the sulcus surrounding the tooth, providing a clean and dry operating field during restorative and prosthodontic procedures.

COMPOSITION

Aluminium Chloride – 25%, Glycerine, Fillers, Colouring Pigment, Flavouring Agents, Purified Water

INTENDED USE

Preva Block is intended for use as a topical hemostatic agent to control minor gingival bleeding and gingival crevicular fluid during dental procedures.

INDICATIONS FOR USE

Preva Block is indicated for the control of minor gingival bleeding and fluids during:

- Restorative procedures.
- Crown and bridge preparations.
- Impression-taking procedures.
- Gingival sulcus management and retraction.
- Oral prophylaxis and other dental procedures where temporary control of minor gingival bleeding is required.

WARNINGS

- For professional dental use only.
- The product should be used only by, or under the supervision of, a qualified dental professional.
- Do not inject.
- For topical dental use only.
- Use only for temporary control of minor gingival bleeding and gingival crevicular fluid during dental procedures.
- Do not swallow.
- Avoid unnecessary or prolonged contact with oral soft tissues.
- Keep out of reach of children.

PRECAUTIONS

- Avoid contact with eyes and skin.
- Do not ingest or take internally.
- Avoid excessive application of the material.
- Rinse the treatment area thoroughly after use.
- Ensure that residual material is removed before proceeding with adhesive or restorative procedures.

Eye Contact

- The product may cause eye irritation. In case of accidental eye contact, immediately rinse thoroughly with plenty of clean water. Seek medical attention if irritation persists or if clinically indicated.

Skin Contact

- Avoid contact with skin. The product may cause irritation or sensitization in susceptible individuals. In case of accidental skin contact, immediately remove the material and wash the affected area thoroughly with soap and water. Seek medical attention if irritation or sensitization develops or persists.

Accidental Ingestion

- If accidentally swallowed, rinse the mouth with water. Do not induce vomiting. Seek medical advice if a significant quantity has been swallowed or if adverse symptoms occur.

CONTRAINDICATIONS

Preva Block should not be used:

- In patients with a known allergy or hypersensitivity to Aluminium Chloride or any other component of the product.
- For the control of bleeding associated with surgical extraction of teeth.
- On areas where use of a topical gingival hemostatic agent is clinically inappropriate.

INSTRUCTIONS FOR USE

1. Preparation

- Remove the protective cap from the syringe.
- Securely attach a new applicator tip to the syringe.
- Ensure that the treatment area is accessible and adequately isolated.

2. Direct Application

- Apply a sufficient amount of Preva Block directly into the gingival sulcus or onto the area of minor gingival bleeding.
- Carefully follow the contour of the prepared tooth and ensure adequate contact with the affected gingival tissue.
- Allow the gel to remain in place for approximately 2 minutes.
- Rinse the treated area thoroughly with water.
- Gently dry and inspect the area for adequate hemostasis.
- Once bleeding and fluid are adequately controlled, proceed with the intended dental procedure or impression.

3. If Bleeding Persists

- If adequate hemostasis has not been achieved:
- Thoroughly rinse away the previously applied gel.
- Reapply Preva Block to the affected area.
- Allow it to remain in place for an additional 2 minutes.
- Rinse thoroughly with water and gently dry.
- Confirm adequate hemostasis before proceeding with the dental procedure.

4. Use with Retraction Cord

For gingival sulcus retraction:

- Apply Preva Block to an appropriate gingival retraction cord or use it in conjunction with the cord according to the selected clinical technique.
- Carefully place the impregnated retraction cord step-by-step into the gingival sulcus using a suitable cord-packing instrument.
- Avoid excessive pressure or trauma to the gingival tissue.
- After the required clinical time, carefully remove the retraction cord.
- Thoroughly rinse the treatment area with water and gently dry before proceeding with the impression or restorative procedure.

Important Application Information

- Use a new applicator tip for each patient.
- Do not return dispensed material to the syringe or original container.
- Material that has already been dispensed or applied must not be reused.
- Remove residual material thoroughly by rinsing before proceeding with bonding, impression-taking, or other subsequent procedures.

STORAGE

- Store at 10°C to 24°C (50°F to 75°F).
- Store in a cool, dry place.
- Protect from excessive heat and direct sunlight.
- Keep away from moisture.
- Close the syringe cap securely immediately after use.
- Do not use after the expiry date stated on the packaging.
- When properly closed and stored under the recommended conditions, opening the syringe does not affect the stated shelf life.

PRESENTATION

Ref. No.: PL-0153 – Preva Block – 25% Aluminium Chloride Hemostatic Gel- 3g Syringe

SHELF LIFE

3 years from the date of manufacture, when stored under the recommended storage conditions.

DISPOSAL

Dispose of used applicator tips, residual material, empty syringes, and packaging in accordance with applicable local regulations and the prevailing waste-disposal requirements of the country in which the product is used.

CONTACT

Prevalent Healthcare LLP, Poongode Road, Karutheni, Malappuram, Kerala – 679339, India.
Ph.: +91 4931257954 Email: info@prevalentonline.com
www.prevalentonline.com

Failure to comply with the specified storage conditions may alter the working characteristics of the material and reduce its shelf life.

The manufacturer is not responsible for any loss of quality resulting from non-compliance with the specified conditions for transportation, storage, and use of this product. The user assumes full responsibility for using the material for any purposes other than those specified by the manufacturer. This material is intended solely for professional dental use. Application must be carried out strictly according to the instructions for use. The manufacturer accepts no liability for damages resulting from failure to follow the provided instructions or for use beyond the stipulated applications. The user is responsible for assessing the material's suitability for any purpose not explicitly stated in this instruction sheet. The descriptions and data provided do not constitute a warranty of attributes and are not legally binding.

This product is specially formulated for use in dentistry and to be used by qualified dental personnel

Symbol for "TEMPERATURE LIMITATION"

Symbol for "KEEP OUT OF THE REACH OF CHILDREN"

Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"

Symbol for "KEEP AWAY FROM SUNLIGHT"

Symbol for "MANUFACTURER NAME & ADDRESS"

Symbol for "CONSULT INSTRUCTION FOR USE"

Symbol for "CATALOGUE NUMBER"

Symbol for "MEDICAL DEVICE" Symbol for "KEEP DRY"

Symbol for "NON STERILE" Symbol for "USE BY"

Symbol for "BATCH CODE" Symbol for "CAUTION"

Document No	PL-IFU-038
Revision	02
Date	16-02-2026